Distributed by: • Distribué par : • Vertrieb: • Distribuido por: • Gedistribueerd door: • Distribuito da: • Distribuert av: • Jälleenmyyjä: • Distribuera av: • Distribuert af: •
 Distribuído por: • Distributer: • Distributor: • Forgalmazza: • Dystrybutor: • Dağıtımci: • Distribútor: • Διανέμεται από: • Дистрибутор: • 分销商: • Distribuido por: • تم التوزيع بواسطة: •
 מופץ על-ידי:

Medtronic

Medtronic MiniMed

18000 Devonshire Street
 Northridge, CA 91325
 USA
 1800 646 4633
 +1 818 576 5555
 www.medtronicdiabetes.com

Medtronic B.V.
 Earl Bakkenstraat 10
 6422 PJ Heerlen
 The Netherlands

(01)00763000246044(10)NA

M968733A23B2_A

M968733A23B2_A

2019-02-07

MiniMed™ Mio™ Advance

Infusion Set

Cathéter

Infusionsset

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Infusieset

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Slangesett

Infusiosetti

Infusionsset

Infusionsæt

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Infuzijski set

Infuzní set

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Infúzyon Seti

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Indications for use

The MiniMed™ Mio™ Advance infusion set is indicated for subcutaneous infusion of medication administered by an external pump.

Description

The MiniMed Mio Advance infusion set has a 90-degree soft cannula. It is delivered ready to use in a pre-loaded insertion device with automatic needle retraction. It is provided sterile and nonpyrogenic.

Intended use

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes mellitus. The infusion set has been tested to be compatible with insulin approved for subcutaneous infusion.

Contraindications

This infusion set is indicated for subcutaneous use only. Do not use the infusion set for intravenous infusion. Do not use the infusion set with blood or blood products.

General instructions

- You can wear the infusion set for a maximum of 72 hours, or per your healthcare professional's instructions.
- If you are using the infusion set for the first time, do the first set-up in the presence of your healthcare professional.

Warnings

- Do not use if the package is opened or damaged. An opened or damaged package may

be contaminated. This can cause an infection. This is a single use product.

- Do not use if the disconnect cover is removed. Use a new infusion set instead. One purpose of the disconnect cover is as a safety feature against accidental firing of the insertion device which could cause injury.
- Do not change the infusion set just before bedtime unless you can check your blood glucose one to three hours afterwards.
- Do not put alcohol, disinfectants, perfumes, deodorants, cosmetics, or other substances with solvents on the infusion set. These substances may damage the infusion set. High or low blood glucose can result if the infusion set is damaged.
- Never point the insertion device toward any body part where insertion is not desired.
- Always rotate your insertion site when changing the infusion set. Reusing an insertion site too often may cause scarring and inaccurate insulin delivery. Review your pump user guide on how to rotate insertion sites.
- Check the insertion site often through the clear window. Improper insertion and maintenance of the insertion site can cause inaccurate insulin delivery, infection, or site irritation. Replace the infusion set at a new site if the soft cannula is not properly inserted.
- Check the cannula housing and tubing for blood. Blood could cause insufficient insulin delivery. This may result in high blood glucose. Replace the infusion set at a new site if there is blood.
- Fill tubing completely with insulin before insertion. Do not leave air in the tubing. An

uncontrolled amount of insulin may be delivered if air is present.

- Check the tubing for clogs or leaks if your blood glucose reading is high. Clogs or leaks may restrict insulin delivery and result in high blood glucose. Replace the infusion set even if you are not sure if there is a clog or leak.
- Do not attempt to clear air or an occlusion in the tubing while it is connected to your body. An uncontrolled amount of insulin may be delivered. This can cause high or low blood glucose. Disconnect the tubing before adjusting it.
- Reuse of the infusion set may cause damage to the cannula or needle and lead to infection, site irritation, or inaccurate medication delivery.
- **For the MiniMed tubing connector:** If insulin or any liquid gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. This may result in the delivery of too little or too much insulin, which can cause hyperglycemia or hypoglycemia. If this occurs, start over with a new reservoir and infusion set.
- For No Delivery Alarms, refer to your pump user guide.

Precautions

- Do not change to a product you have not used before without consulting a healthcare professional for correct handling. The infusion set product is closely related to insertion location. Therefore, your healthcare professional should always be included when choosing an infusion set product.
- An inaccurate insertion site or inappropriate site management can result in physical pain,

influence insulin absorption or cause damage to the product. Other site issues may also occur.

- Use an insertion site as recommended by your healthcare professional. The choice of site depends on treatment and patient-specific factors. Check that the insertion site is free of skin irritation such as redness, scar tissue, bleeding, and so forth. Do not insert into muscle or over bone. Doing so can cause pain or damage to the product. Replace the infusion set and use a new insertion site if this happens.
- If needed, remove body hair around the insertion site to ensure that the adhesive sticks to the skin.
- To avoid stopping insulin delivery that could result in high blood glucose, do not use an insertion site where clothing can push or pull on the infusion set. An insertion site such as under the waistband of your pants may cause various set problems including (but not limited to) disconnecting the tubing, pulling out the cannula, or kinking the tubing, any of which could interrupt insulin delivery.
- Do not reposition the infusion set on your body. Repositioning the infusion set may damage the adhesive. Replace the infusion set if the adhesive is damaged.
- Dispose of the insertion device in a proper sharps container according to local laws. The insertion device has a needle. Improper disposal may cause an accidental needle stick.
- Avoid touching a used soft cannula. A used infusion set is infectious waste. Dispose of a used infusion set in a proper container according to local laws.

- Check your blood glucose while disconnected from the pump. Consult your healthcare professional on how to correct for missed insulin and for how much time the pump can remain disconnected.
- Always check your blood glucose one to three hours after inserting a new infusion set. This is to confirm accurate insulin delivery. High or low blood glucose can result if insulin delivery is not accurate.
- Check your blood glucose several times throughout the day or as recommended by your healthcare professional.
- Carefully monitor your blood glucose level when disconnected and after reconnecting your infusion set.
- Store and handle insulin according to the manufacturer's instructions.

Storage

Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight. Do not store infusion sets where it is humid.

Additional information for the healthcare professional

- Consider treatment and patient-specific factors when choosing an insertion site.
- The length of the soft cannula should be based on treatment and patient factors such as physiology and activity level. If the soft cannula is too long, insertion pain or insertion into the bone or muscle may occur. If the soft cannula is too short, leakage or irritation at the insertion site may occur. In both cases, the infusion must

be stopped and a new soft cannula length has to be inserted at a new insertion site.

- Inform the patient about the correct choice of insertion site in relation to the choice of soft cannula length.


Warranty

For product warranty information, please contact your local help line, representative, or visit:

www.medtronicdiabetes.com/warranty

©2019 Medtronic MiniMed, Inc. All rights reserved. MiniMed and Mio are trademarks of Medtronic MiniMed, Inc.

Components

- | | |
|--|--|
| A. Top button | H. Paper backing |
| B. Disconnect cover | I. Adhesive |
| C. Insertion device | J. MiniMed or Luer-lock tubing connector |
| D. Raised arrows indicating tubing direction | K. Tubing |
| E. Clear Window | L. Site Connector |
| F. Cannula housing | M. White cap for site connector |
| G. Soft cannula | |

Instructions for use

- Read all instructions carefully before using the MiniMed Mio Advance infusion set.
- Read your pump user guide before connecting the infusion set to your pump for information about pump therapy. This includes information

on connection and priming procedures, possible errors, and potential risks and their causes.

- You should follow good hygienic procedures. If you are using this infusion set for the first time, do the first set-up in the presence of your healthcare professional.
- Do not insert the product at temperatures below 10°C (50°F).

Inserting

Unfold the front cover of this booklet to view images that are keyed to the steps below.

1. Wash your hands with soap and water.
2. Select an insertion site on a recommended area (shown in gray) and as recommended by your healthcare professional.
3. Clean the insertion site with a disinfectant as directed by your healthcare professional. Allow to air dry before inserting the infusion set. Remove hair around the insertion site to ensure the adhesive sticks to the skin if needed.
4. Open the package. Then, remove the paper from the tubing.
5. Place the tubing connector on the top of a newly filled reservoir. Twist the MiniMed tubing connector clockwise until it locks in place or twist the Luer-lock tubing connector until it is tightly secured.

WARNING: For the MiniMed tubing connector, make sure there is no liquid inside the tubing connector or on top of the reservoir. Liquid can block vents and cause inaccurate insulin flow. This can

result in high or low blood glucose. If there is liquid, use a new infusion set and a new reservoir.

6. Remove the white cap from the site connector by (a) squeezing the sides of the site connector then (b) pulling the white cap off. Keep the white cap for later use.
7. Put the reservoir into your pump, then completely fill the tubing with insulin as instructed by your pump user guide. Do not leave any air in the tubing. The tubing is completely filled when you see drops of insulin coming out of the site connector.
8. Remove the paper backing from the adhesive. Be careful not to touch the adhesive.
9. Remove the disconnect cover from the insertion device by (a) gently squeezing the sides of the disconnect cover then (b) pulling it away from the insertion device. Keep the disconnect cover for later use. The raised arrows indicate the tubing direction.
- 10a. Stretch the skin until smooth. Then press the insertion device against your skin.
- 10b. Press the top button completely down to insert the infusion set.
11. Gently and carefully remove the insertion device.
12. Press the adhesive onto the skin with your finger. Replace the infusion set if the adhesive does not stick to the skin.
13. Gently hold the cannula housing steady with your finger. Then push the site connector

straight into the cannula housing until you hear a click.

14. Fill the soft cannula with insulin:

6 mm: 0.6 units (0.006 ml)

9 mm: 0.6 units (0.006 ml)

Dispose of the insertion device in an appropriate sharps container according to local laws.

Keep the disconnect cover and the white cap for use when disconnecting the infusion set.

Disconnecting

- Wash your hands with soap and water.
- Gently hold the cannula housing steady with your finger. Then (a) squeeze the sides of the site connector and (b) pull it out from the cannula housing.
- Put the white cap on the site connector.
- Put the disconnect cover on the cannula housing. Then push the cover into the cannula housing until you hear a click.

Reconnecting

- Wash your hands with soap and water.
- Use your finger to gently hold the cannula housing steady. Then (a) squeeze the sides of the disconnect cover and (b) pull it out from the cannula housing.
- Remove the white cap from the site connector by (a) squeezing the sides of the site connector then (b) pulling the white cap off. Ensure that there is no air in the tubing. **ONLY if there is air in the tubing:** Fill the tubing with insulin as instructed in your pump user

guide. Do not leave any air in the tubing.

- Use your finger to gently hold the cannula housing steady. Then, push the site connector into the cannula housing until you hear a click.

Removing

- Wash your hands with soap and water.
- Use your finger to gently hold the cannula housing steady. Then (a) squeeze the sides of the site connector and (b) pull it out from the cannula housing.
- Disconnect the tubing from the pump and remove the tubing connector from the reservoir by turning the connector counter-clockwise. For specific instructions to disconnect tubing from the pump, please refer to your pump user guide.
- Carefully lift the adhesive around the cannula housing. Then, pull the soft cannula out of the skin.

Dispose of your used infusion set in an appropriate container according to local laws.

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Medtronic (Taiwan) Ltd.
Medtronic (Thailand) Ltd.
Ticaret Ltd. Sirketi.
Medtronic Ltd.
Medtronic Österreich GmbH