
One-press Serter

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Introduction

The one-press serter (MMT-7512) is used as an aid for inserting Medtronic MiniMed® glucose sensors.

- A. Bump on both buttons
- B. Thumbprint marking

Potential risks related to serter use

General risks with serter use may include skin infection around the area where the serter is used.

Indications for use

The serter is used as an aid for inserting the sensor. It is indicated for single-patient use and is not intended for multiple patient use.

Contraindications

Do not use the serter on products other than the approved Medtronic MiniMed glucose sensors. Medtronic cannot guarantee the safety or efficacy of this product if used with other products.

General warnings

Read this entire user guide before attempting to insert the sensor. The serter does not work the same as other Medtronic insertion devices. Failure to follow directions may result in improper insertion, pain, or injury.

Contact the 24 Hour HelpLine or your local representative, if you have questions or concerns.

Do not make therapy decisions based on sensor glucose values because sensor glucose and blood glucose values may differ. If your sensor glucose reading is low or high, or if you feel symptoms of low or high glucose, do the following prior to making therapy decisions. Confirm your blood glucose with your BG meter using a fingerstick blood sample.

A retractable needle is attached to the sensor and minimal blood splatter may occur. If you are a healthcare professional or caregiver, wrap sterile gauze around the sensor to minimize contact with blood. Keep as much distance as possible between you and the patient when removing the needle.

Keep the needle housing within sight at all times to avoid an accidental needlestick or puncture.

Always inspect the sensor packaging for damage before use. Sensors are sterile and non-pyrogenic, unless the package has been opened or damaged. Do not use the sensor if the sterile package has been opened or damaged. Use of an unsterile sensor can cause site infection.

Watch for bleeding at the insertion site (under, around, or on top of the sensor).

If bleeding occurs, do the following:

1. Apply steady pressure using sterile gauze or a clean cloth, placed on top of the sensor for up to three minutes. The use of unsterile gauze can cause site infection.

2. If bleeding stops, connect the transmitter to the sensor.

If bleeding does not stop, do not connect the transmitter to the sensor because blood can get into the transmitter connector, and could damage the device.

If bleeding continues, do the following:

1. Remove the sensor and continue to apply steady pressure until the bleeding stops. Discard the sensor in a sharps container.
2. Check the site for redness, bleeding, irritation, pain, tenderness, or inflammation. Treat based on instructions from your healthcare professional.
3. Insert a new sensor in a different location.

Make sure the sensor is securely placed in theserter to avoid improper insertion, pain, or minor injury. Keep theserter away from children. This product contains small parts and may pose a choking hazard. Refer to the sensor user guide for additional sensor related warnings.

General precautions

Wash your hands with soap and water before inserting the sensor to help prevent site infection.

If you are a healthcare professional or caregiver, wear gloves when inserting the sensor into someone other than yourself to avoid contact with patient blood.

Do not insert the sensor through tape. Inserting the sensor through tape may cause improper sensor insertion and function.

Only use alcohol to prepare the insertion site, to ensure that residue is not left on the skin.

Rotate the sensor insertion sites so that they do not become overused.

Discard used sensors and needle housings in a sharps container after each use to avoid accidental needlestick or puncture.

Do not clean, resterilize, or try to extract the needle from the needle housing. An accidental needlestick or puncture may occur.

Do not reuse sensors. Reuse of a sensor may cause damage to the sensor surface and lead to inaccurate glucose values, site irritation, or infection.

Always make sure the sensor is firmly attached so that it does not come out of your body. Factors that may increase the likelihood of the sensor falling off include increased physical activity (particularly for younger patients) and improper taping technique.

Where to insert the sensor

CAUTION: Avoid the 2 inch (5.0 cm) area around the navel to help ensure a comfortable insertion site and to help with sensor adhesion.

Choose an insertion site that has an adequate amount of subcutaneous fat. For a diagram of the sensor insertion area, see the sensor user guide.

Do not insert the sensor in muscle or areas constrained by clothing or accessories, areas with tough skin or scar tissue, sites subjected to rigorous movement during exercise, or in sites under a belt or on the waistline for best sensor performance and to avoid accidental sensor removal.

Inserting the sensor

Always refer to the sensor user guide for instructions on how to insert the sensor.

Storage and handling

Store the serter in the released position, at room temperature, to maintain optimum product performance and life.

Specifications

Approximate Dimensions	Approximate Weight
3.09 x 2.72 x 2.27 inches (7.85 x 6.91 x 5.77 centimeters)	3.97 ounces (112.5 grams)

Serter life of use

The serter has a minimum service life of one year.

Serter maintenance

The serter is intended for single-patient use. The serter needs to be cleaned after every use or whenever there is debris, blood, or other contaminants.


Cleaning

To clean the serter, you will need the following materials: mild liquid soap (for example, Ivory® liquid soap), a soft-bristled toddler toothbrush, and a container.

To clean the serter:

1. Prepare a mild liquid soap solution using 1 teaspoon (5 milliliters) of mild liquid soap per 1 gallon, which is equal to 16 cups (3.8 liters) of tap water.
2. Rinse the serter under running tap water at room temperature for at least one minute. Continue rinsing until the serter is visibly clean. Make sure all hard-to-reach areas are rinsed completely.
3. When rinsing the serter, press and release the mechanism on the underside of the serter to ensure that the entire device is completely rinsed with water.
4. Submerge the serter in the mild liquid soap solution and soak it for at least 10 minutes.
5. Brush the entire surface of the serter using a soft-bristled toddler toothbrush until visibly clean. Pay close attention to brushing the hard-to-reach areas.
6. Rinse the serter under running tap water at room temperature for at least one minute, and until all visible liquid soap is gone.
7. Shake off any excess water or moisture and place the serter upright on a clean, dry, cloth to air dry.

Icon table

	Consult instructions for use.		Caution: Specific warnings may not be found on label.
	Batch code		Humidity limitation
	One per container/package		Date of manufacture (YYYY-MM-DD)
	Temperature limit		Prescription only
	Catalogue or model number		Manufacturer

©2015, Medtronic MiniMed, Inc. All rights reserved. MiniMed® is a trademark of Medtronic MiniMed, Inc.

Ivory® is a registered trademark of The Procter & Gamble Company.

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